

## **Blake**

## **Cushion Comparison Chart**

			GEL	GEL			FOAM & GEL							
_	BRONZE	SILVER	GOLD	GOLD+	BOLD	PC Cutback	VINGA	HYBRID	HYBRID LP	HYBRID DEEP	HY. DEEP LP	SPLIT	PLATINUM	PLAT. PLUS
Great For:	*	* () N	* ~ ~	* ~ ~	* ~ ~	()	141	† <b>‡</b> ↓ <b>▽</b>	Ü	世へ	() N	144 (7 🚓	* 🕏	* ~
Risk of Skin Breakdown	*	***	***	***	***	*	**	***	**	***	**	***	***	***
Shape/ Contour	*	*	*	*	*	**	**	**	**	***	***	**	**	**
Immersion	**	***	***	***	***	***	***	**	**	***	**	**	***	***
Envelopment	**	***	***	***	***	**	***	**	**	***	***	**	***	***
Stability	**	**	**	**	**	**	**	***	***	***	***	***	**	**
Shear reduction	***	***	***	***	***	*	**	**	**	**	**	**	***	***
Vibration Dampening	**	***	***	***	***	*	*	**	**	**	**	**	***	***
Adjustment	***	***	***	***	***		*	**	**	**	**	***	**	**
Temp Control	**	**	**	**	**	***	***	*	*	*	*	*	*	*
Product Weight	**	***	***	***	***	*	*	**	**	**	**	**	***	***
Weight Capacity	*	*	*	**	***	*	*	*	*	*	*	*	**	***

























Blake Medical realizes that everyone is different and requires different cushion attributes to make each individual comfortable and safe. Please note that the above chart is for reference only and that each individual user should contact a healthcare professional to ensure which cushion is right for them.